

XPANSION[®]

MICRO-AUTOGRAFTING KIT

Dressing Protocol



Post-Operative Care

- Non-adherent perforated dressing can be used (e.g. bridal veil) as contact layer.
- Secure non-adherent dressing to wound with Steri-Strips.
- Maintain a moist environment (e.g. hydrogel or moisture-holding dressing).
- Splint or immobilize if in mobile area.
- Outer dressing can be changed at 5 - 7 days leaving contact layer intact.
- Irrigation and/or soak can be used at dressing changes.
- Contact dressing can be gently teased off graft at 10 - 14 days using saline irrigation to prevent any dressing adherence.



Applied Tissue
Technologies

99 Derby St #200
Hingham, MA 02043
www.appliedtissue.com
(781) 366-3848